

# Antibiotic purification

Adsorbent and ion exchange resins can be used in the purification process of antibiotics. They are well established in the purification of molecules like cephalosporins, glycopeptides and aminoglycosides.

## Cephalosporins


Cephalosporins are the most important antibiotics having  $\beta$ -lactam ring. They are obtained from the fungus *Acremonium chrysogenum*, also known as *cephalosporium* and they have a wide use against bacteria in various severe infections such as respiratory tract infection (RTI), skin infection and urinary tract infection (UTI). 7 Aminocephalosporanic acid (7-ACA) is made from Cephalosporin C (CPC) and is a key intermediate for synthesizing the four major classes of cephalosporin antibiotics.

After fermentation, the separation of biomass and antibiotic-containing broth is generally achieved by microfiltration by which the biomass is removed from the CPC containing filtrate. The filtered broth is then passed through large-scale hydropho-

bic interaction chro-

matography (HIC) columns to remove impurities, in particular proteins, peptides, salts and side products like deacetyl CPC (DAC) and deacetoxy CPC (DAOC).


The first column, called a scavenger, is filled with an adsorbent, e.g. **AMBERLITE™ XAD™ 1180N**. The second column is filled with another adsorbent, e.g. **AMBERLITE™ XAD™ 1600N**. The remaining color is removed by percolating the CPC solution through a column filled with a weak anionic acrylic resin in acetate form, **AMBERLITE™ FPA53**. Then CPC is converted into 7-ACA by a two step cleavage with D-amino acid oxidase (DAO) and glutaryl acylase (GAC). A Merry Go Round column system is typically used for this process.



## Vancomycin

Vancomycin is one of the glycopeptide antibiotics. It has bactericidal activity against aerobic and anaerobic Gram-positive bacteria and is used in the treatment of methicillin resistant staphylococcus aureus (MRSA). The initial purity of vancomycin in a fermentation broth of bacteria *S. orientalis* is around 32–35%, and the purity must be driven up to greater than 95% for pharmaceutical application.

The capture step can be performed on a strongly acidic resin, **AMBERLITE™ FPC23H**, decolorization and desalting of crude filtered vancomycin can be performed on **AMBERLITE™ XAD™16** or **XAD™1600N** and the final purification step on **AMBERCHROM™ CG161M** enables >95% product purity.



## Tobramycin

Tobramycin is an aminoglycoside antibiotic derived from *Streptomyces tenebrarius*. It is used to treat various types of bacterial infections, particularly Gram-negative infections. The product of microbial fermentation is carbamoyltobramycin, which is converted to tobramycin using ammonium hydroxide hydrolysis. Tobramycin is then recovered upon a carboxylic resin such as **AMBERLITE™ FPC3500** in the ammonium form and then decolorized upon an anionic resin such as **AMBERLITE™ FPA90 CL** or **AMBERLITE™ FPA40 CL**. A further separation of A, B and C components can be performed upon **AMBERLITE™ CG50 TYPE 1**. A similar process exists to produce other aminoglycoside antibiotics like Kanamycin, Netilmicin, Sisomicin or Gentamicin.

Biomolecule	Application	Process steps	Resin
Cephalosporins	β-lactam for respiratory tract infection, Skin Infection, Urinary Infection	Scavenger Extraction Decolorization	AMBERLITE™ XAD™ 1180N AMBERLITE™ XAD™ 1600N AMBERLITE™ FPA53
Vancomycin	Glycopeptide for Gram-positive bacteria MRSA	Capture Decolorization & desalting Final purification	AMBERLITE™ FPC23H AMBERLITE™ XAD™16N AMBERCHROM™ CG161M
Tobramycin	Aminoglycosides for gram negative infections	Capture Decolorization Final purification	AMBERLITE™ FPC3500 AMBERLITE™ FPA90 CL AMBERLITE™ CG50 (Type 1)
Streptomycin	Aminoglycosides for treatment of active tuberculosis	Conversion to sulfate form	AMBERLITE™ FPA40
Oritavancin	Glycopeptide for Gram-positive bacteria	Capture Decolorization & desalting Final purification	AMBERLITE™ FPC23H AMBERLITE™ XAD™ 16N AMBERCHROM™ CG161M
Erythromycin	Macrolide antibiotic for respiratory tract infection, Skin Infection, STD	Capture and concentration Decolorization Desalting	AMBERLITE™ FPC3500 AMBERLITE™ FPA98 CL AMBERLITE™ XAD™ 16N
Geldanamycin	Ansamycin for broad spectrum antibiotic	Capture	AMBERLITE™ XAD™ 1600N

Picture credits p.1: istock, p.2: DuPont Europe GmbH



Water Solutions

[dupont.com/water](http://dupont.com/water)

NA: 1 800 447 4369

DuPont™, the DuPont Oval Logo, AMBERCHROM™, AMBERLITE™ are trademarks or registered trademarks of DuPont or its affiliates. Copyright © 2019 E. I. du Pont de Nemours and Company.